

Teklad Diamond Twists

Teklad 7979C, 7979HC, 7978C



Certified Irradiated Diamond Twists now come in 3 different sizes

7979C – Full sized Certified Irradiated Diamond Twists, approximately 10 inches in length

7979HC – Diamond Half Twists are approximately 5 inches in length

NEW!!

7978C – Diamond Twist Minis approximately 2 inch pieces in length. Excellent foraging material!



All of our Twists can be used as both enrichment and a source of nesting material. Our new and unique Diamond Twist Minis can also be used as foraging material, as well as enrichment and nesting.

Teklad 7978C Diamond Twists Minis are the first environmental and foraging material that is both Certified and Irradiated.

PRODUCT DESCRIPTION

Made from thick rolled tissue paper with very low dust content, Diamond Twists can be used conveniently in GLP toxicology studies and are ideal for use with nude mice. In multiple studies this product experienced 100% acceptance for shredding and did not wind around rodents' legs or pups. The paper material will not clog drains or disposals.

PRODUCT FEATURES AND BENEFITS

- First and only environmental enrichment material that is both certified and irradiated
- Can be autoclaved in the cage
- Highly defined raw material
- Monitored for quality and consistency throughout production

- Produced in a plant dedicated to the production of lab animal bedding
- Premium white color aids in observation for health monitoring
- 100% acceptance and for shredding
- Does not wind around rodents' legs, or pups
- Won't clog drains or disposals
- Product come in 3 sizes
- **NEW** Diamond Twist Minis are sized to use as foraging material

PRODUCT FORMATS

Certified Irradiated Diamond Twists are monitored for quality and consistency throughout production, and the material is routinely certified for pesticides and heavy metals. It is a certified, irradiated enrichment product with certificates of analysis available on a per-lot basis.

PACKAGING

All Diamond Twists products are certified and irradiated. Product can be autoclaved in cages or other packaging.

Product Name	Product Number	Packaging
Diamond Twists	7979C	1900 pieces/case Case contains 9 x 2 lb vacuum packages
Diamond Half Twists	7979HC	3700 pieces/case Case contains 9 x 2lb vacuum packages.
Diamond Twist Minis	7978C	22 lb bag (10 kg) poly lined sealed bag

Find out more Contact your Inotiv representative or visit us at inotiv.com