

# TD.00502

## Rodent Diet (2018, 200 Doxycycline)

Formula	g/Kg
2018, Teklad Global 18% Protein Rodent Diet	999.8
Doxycycline Hyclate	0.2

### Footnote

This diet was designed to deliver a daily dose of 0.5-1 mg of doxycycline based on consumption of 3-5 g/d by a mouse. Doxycycline hyclate contains approximately 87% doxycycline.

### Selected Nutrient Information <sup>1</sup>

	% by weight	% kcal from
Protein	18.6	24.4
Carbohydrate <sup>2</sup>	44.2	57.9
Fat	6.0	17.7
<b>Kcal/g</b>	<b>3.1</b>	

<sup>1</sup> Values are calculated from ingredient analysis or manufacturer data

<sup>2</sup> Estimated digestible carbohydrate

### Collaborate with a Nutritionist

- [askanutritionist@inotivco.com](mailto:askanutritionist@inotivco.com)

Teklad diets are designed & manufactured for research purposes only

Teklad™ Custom Diets, Madison, WI | [inotivco.com](http://inotivco.com)



### Key Features

- Standard Diet Base
- Doxycycline
- Gene Induction
- 2018

### Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week  
*More frequent replacement may be advised*
- Lead time:
  - 2 weeks non-irradiated
  - 4 weeks irradiated



### Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

### Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

### Place Your Order (USA & Canada)

Please Choose One

- [www.inotivco.com/teklad-orders](http://www.inotivco.com/teklad-orders)
- [tekladorders@inotivco.com](mailto:tekladorders@inotivco.com)
- 800.483.5523

### Price Quotes | Order Status

- [teklad@inotivco.com](mailto:teklad@inotivco.com)