

# TD.120528

42% Kcal/Fat Diet (Incr. Sucrose, 1.25% Chol.)

Formula	g/Kg
Casein	195.0
DL-Methionine	3.0
Sucrose	405.36
Maltodextrin	75.0
Anhydrous Milkfat	210.0
Cholesterol	12.5
Cellulose	50.0
Mineral Mix, AIN-76 (170915)	35.0
Calcium Carbonate	4.0
Vitamin Mix, Teklad (40060)	10.0
Ethoxyquin, antioxidant	0.04
Green Food Color	0.1

## Footnote

Modified from TD.10364 to increase added cholesterol to 1.25% (total ~1.3%). Originally modified from TD.88137 for increased sucrose.

## Selected Nutrient Information <sup>1</sup>

	% by weight	% kcal from
Protein	17.3	15.2
Carbohydrate	48.5	42.7
Fat	21.2	42.0
<b>Kcal/g</b>	<b>4.5</b>	

<sup>1</sup> Values are calculated from ingredient analysis or manufacturer data

## Collaborate with a Nutritionist

- [askanutritionist@inotivco.com](mailto:askanutritionist@inotivco.com)

Teklad diets are designed & manufactured for research purposes only

Teklad™ Custom Diets, Madison, WI | [inotivco.com](http://inotivco.com)



## Key Features

- Purified Diet
- Modification of TD.88137
- Increased Sucrose
- High Cholesterol

## Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week  
*More frequent replacement may be advised*
- Lead time:
  - 2 weeks non-irradiated
  - 4 weeks irradiated



## Product Specific Information

- 1/2" Pellet or Powder (crumbly)
- Minimum order 3 Kg
- Irradiation available upon request

## Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

## Place Your Order (USA & Canada)

Please Choose One

- [www.inotivco.com/teklad-orders](http://www.inotivco.com/teklad-orders)
- [tekladorders@inotivco.com](mailto:tekladorders@inotivco.com)
- 800.483.5523

## Price Quotes | Order Status

- [teklad@inotivco.com](mailto:teklad@inotivco.com)