

# TD.140263

Doxycycline Diet (625, 2016, Blue)

## Formula

	g/Kg
2016, Teklad Global 16% Protein Rodent Diet	999.075
Doxycycline Hyclate	0.625
Blue Food Color	0.3

## Footnote

Diet adds 625 mg/kg doxycycline hyclate to Teklad's Global 2016 Rodent Diet and is color coded blue. Diet was designed to deliver a daily dose of 2 - 3 mg of doxycycline based on consumption of 3 - 5 g/d by a mouse. Doxycycline hyclate contains approximately 87% doxycycline.

## Selected Nutrient Information <sup>1</sup>

	% by weight	% kcal from
Protein	16.2	21.1
Carbohydrate <sup>2</sup>	52.8	68.6
Fat	3.5	10.4
<b>Kcal/g</b>	<b>3.1</b>	

<sup>1</sup> Values are calculated from ingredient analysis or manufacturer data

<sup>2</sup> Estimated digestible carbohydrate

## Collaborate with a Nutritionist

- [askanutritionist@inotivco.com](mailto:askanutritionist@inotivco.com)

Teklad diets are designed & manufactured for research purposes only

Teklad™ Custom Diets, Madison, WI | [inotivco.com](http://inotivco.com)



## Key Features

- Standard Diet Base
- Doxycycline
- Teklad's Globa 2016
- Color Coded Blue

## Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week  
*More frequent replacement may be advised*
- Lead time:
  - 2 weeks non-irradiated
  - 4 weeks irradiated



## Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

## Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

## Place Your Order (USA & Canada)

Please Choose One

- [www.inotivco.com/teklad-orders](http://www.inotivco.com/teklad-orders)
- [tekladorders@inotivco.com](mailto:tekladorders@inotivco.com)
- 800.483.5523

## Price Quotes | Order Status

- [teklad@inotivco.com](mailto:teklad@inotivco.com)