TD.160785

52 kcal/Fat Diet (C16:0, HVO, AMF, Choline/Met)

Formula	g/Kg
Casein, "Vitamin-Free" Test	140.0
Sucrose	341.46
Maltodextrin	119.6156
Hydrogenated Vegetable Oil	190.0
Anhydrous Milkfat	60.0
Palmitic Acid	40.0
Cholesterol	12.5
Cellulose	50.0
Mineral Mix, AIN-76 (170915)	35.0
Calcium Carbonate	4.0
Vitamin Mix, w/o choline, A, D, E (83171)	5.0
Vitamin E, DL-alpha tocopheryl acetate (500 IU/g)	0.1
Vitamin A Palmitate (500,000 IU/g)	0.04
Vitamin D ₃ , cholecalciferol (500,000 IU/g)	0.0044
Choline Dihydrogen Citrate	2.28



Key Features

- Purified Diet
- Rodent NASH Diet
- Cholesterol
- Palmitic Acid

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
 More frequent replacement may be advised
- Lead time:
 - · 2 weeks non-irradiated
 - · 4 weeks irradiated



Product Specific Information

- Powder (crumbly)
- Minimum order 3 Kg
- Irradiation available upon request

Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Selected Nutrient Information

ociceted itatifelit illioilliation		
	% by weight	% kcal from
Protein	12.6	10.2
Carbohydrate	46.2	37.3
Fat	29.0	52.6

A high fat and high sucrose diet with 1.25% added cholesterol. Fat sources include

anhydrous milkfat, Hydrogenated Vegetable Oil and palmitic acid. Hydrogenated vegetable

oil is added as a source of dietary trans fats. No supplemental Met is added. This formula

with VFT casein is no longer available in pellet form. You can order in powder or "dough" form using product code TD.160785.PWD or order TD.190142 for version with regular

Kcal/g 5.0

casein that can be pelleted.

Footnote

Values are calculated from ingredient analysis or manufacturer data

Collaborate with a Nutritionist

askanutritionist@inotivco.com

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotivco.com
- 800.483.5523

Price Quotes | Order Status

teklad@inotivco.com