

TD.180939

40% Kcal Fat (HVO, 2% Chol, 20% Kcal/Fruct)

Formula	g/Kg
Casein	200.0
L-Cystine	3.0
Fructose	225.0
Sucrose, fine ground	125.0
Maltodextrin	130.0
Hydrogenated Vegetable Oil	155.0
Soybean Oil	25.0
Lard	20.0
Cholesterol	20.0
Cellulose	50.0
Mineral Mix, AIN-93G-MX (94046)	35.0
Vitamin Mix, AIN-93-VX (94047)	10.0
Choline Bitartrate	2.0

Footnote

A diet formulated with approximately 40% kcal from fat, 2% cholesterol and 20% kcal from fructose and commonly referred to as the AMLN (Amylin Liver NASH) diet. Hydrogenated vegetable oil is added as a source of dietary trans fats.

Selected Nutrient Information ¹

	% by weight	% kcal from
Protein	17.7	15.8
Carbohydrate	49.1	43.7
Fat	20.2	40.5
Kcal/g	4.5	

¹ Values are calculated from ingredient analysis or manufacturer data

Collaborate with a Nutritionist

- askanutritionist@inotivco.com

Teklad diets are designed & manufactured for research purposes only
Teklad™ Custom Diets, Madison, WI | inotivco.com



Key Features

- Purified Diet
- Hydrogenated Vegetable Oil
- Cholesterol
- Fructose

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (crumbly)
- Minimum order 3 Kg
- Irradiation not advised
Contact a nutritionist for recommendations

Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotivco.com
- 800.483.5523

Price Quotes | Order Status

- teklad@inotivco.com