

# TD.190142

52 kcal/Fat Diet (C16:0, HVO, AMF, Choline/Met)

Formula	g/Kg
Casein	140.0
Sucrose	341.46
Maltodextrin	119.6156
Hydrogenated Vegetable Oil	190.0
Anhydrous Milkfat	60.0
Palmitic Acid	40.0
Cholesterol	12.5
Cellulose	50.0
Mineral Mix, AIN-76 (170915)	35.0
Calcium Carbonate	4.0
Vitamin Mix, w/o choline, A, D, E (83171)	5.0
Vitamin E, DL-alpha tocopheryl acetate (500 IU/g)	0.1
Vitamin A Palmitate (500,000 IU/g)	0.04
Vitamin D <sub>3</sub> , cholecalciferol (500,000 IU/g)	0.0044
Choline Dihydrogen Citrate	2.28

## Footnote

A high fat and high sucrose diet with 1.25% added cholesterol. Fat sources include anhydrous milkfat, Hydrogenated vegetable oil (source of trans fats) and palmitic acid. No supplemental Met is added. This diet uses the same formula as TD.140154/TD.160785, but uses regular casein instead of "vitamin free test" casein. The change should have no nutritional impact while offering a benefit of better formed pellets.

## Selected Nutrient Information <sup>1</sup>

	% by weight	% kcal from
Protein	12.2	9.8
Carbohydrate	46.2	37.3
Fat	29.1	52.9
Kcal/g	5.0	

<sup>1</sup> Values are calculated from ingredient analysis or manufacturer data



## Collaborate with a Nutritionist

- [askanutritionist@inotivco.com](mailto:askanutritionist@inotivco.com)

Teklad diets are designed & manufactured for research purposes only

Teklad™ Custom Diets, Madison, WI | [inotivco.com](http://inotivco.com)



## Key Features

- Purified Diet
- Rodent NASH Diet
- Cholesterol
- Palmitic Acid

## Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week  
*More frequent replacement may be advised*
- Lead time:
  - 2 weeks non-irradiated
  - 4 weeks irradiated

## Product Specific Information

- 1/2" Pellet or Powder (crumbly)
- Minimum order 3 Kg
- Irradiation available upon request

## Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

## Place Your Order (USA & Canada)

Please Choose One

- [www.inotivco.com/teklad-orders](http://www.inotivco.com/teklad-orders)
- [tekladorders@inotivco.com](mailto:tekladorders@inotivco.com)
- 800.483.5523

## Price Quotes | Order Status

- [teklad@inotivco.com](mailto:teklad@inotivco.com)