

TD.200591

40% Kcal Fat (Palm Oil, 2% Chol, 20% Kcal/Fruct)

Formula	g/Kg
Casein	200.0
L-Cystine	3.0
Fructose	225.0
Sucrose, fine ground	125.0
Maltodextrin	125.0
Cellulose	50.0
Palm Oil	155.0
Soybean Oil	25.0
Lard	20.0
Cholesterol	20.0
Mineral Mix, AIN-93G-MX (94046)	35.0
Vitamin Mix, AIN-93-VX (94047)	15.0
Choline Bitartrate	2.0

Footnote

A diet formulated with approximately 40% kcal from fat, 2% cholesterol and 20% kcal from fructose. Fat sources include palm oil, soybean oil and lard. Designed with similarities to Research Diets, Inc. formula D09100310 which is commonly referred to as the GAN (Gubra Amylin NASH) diet. Vitamins are increased to make the diet more suitable for irradiation.

Selected Nutrient Information ¹

	% by weight	% kcal from
Protein	17.7	15.8
Carbohydrate	49.1	43.7
Fat	20.2	40.5
Kcal/g	4.5	

¹ Values are calculated from ingredient analysis or manufacturer data

Collaborate with a Nutritionist

- askanutritionist@inotivco.com

Teklad diets are designed & manufactured for research purposes only
Teklad™ Custom Diets, Madison, WI | inotivco.com



Key Features

- Purified Diet
- Cholesterol
- Fructose
- Palm Oil

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (crumbly)
- Minimum order 3 Kg
- Irradiation available upon request

Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotivco.com
- 800.483.5523

Price Quotes | Order Status

- teklad@inotivco.com