

# TD.89247

## 60% Fructose Diet

Formula	g/Kg
Casein	207.0
DL-Methionine	3.0
Fructose	600.0
Lard	50.0
Cellulose	79.81
Mineral Mix, Rogers-Harper (170760)	50.0
Zinc Carbonate	0.04
Vitamin Mix, Teklad (40060)	10.0
Green Food Color	0.15

### Footnote

This formula is a modification of diet TD 78463 to increase zinc levels. The vitamin mix contributes 0.4% carbohydrate from the corn starch diluent. Fructose is hygroscopic, and the diet will become tacky after a few days of exposure to the ambient environment.

### Selected Nutrient Information <sup>1</sup>

	% by weight	% kcal from
Protein	18.3	20.2
Carbohydrate	60.4	66.8
Fat	5.2	13.0
Kcal/g	3.6	

<sup>1</sup> Values are calculated from ingredient analysis or manufacturer data

### Collaborate with a Nutritionist

- [askanutritionist@inotivco.com](mailto:askanutritionist@inotivco.com)

Teklad diets are designed & manufactured for research purposes only  
Teklad™ Custom Diets, Madison, WI | [inotivco.com](http://inotivco.com)



### Key Features

- Purified Diet
- Fructose
- Color Coded
- Rodent

### Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week  
*More frequent replacement may be advised*
- Lead time:
  - 2 weeks non-irradiated
  - 4 weeks irradiated



### Product Specific Information

- 1/2" Pellet or Powder (crumbly)
- Minimum order 3 Kg
- Irradiation available upon request

### Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

### Place Your Order (USA & Canada)

Please Choose One

- [www.inotivco.com/teklad-orders](http://www.inotivco.com/teklad-orders)
- [tekladorders@inotivco.com](mailto:tekladorders@inotivco.com)
- 800.483.5523

### Price Quotes | Order Status

- [teklad@inotivco.com](mailto:teklad@inotivco.com)