

TD.96348

Diet (20% Lactose, 2.0% Ca, 1.25% P)

Formula	g/Kg
Casein, "Vitamin-Free" Test	192.0
DL-Methionine	3.0
Sucrose	264.78
Lactose, monohydrate	200.0
Corn Starch	150.0
Corn Oil	50.0
Cellulose	50.0
Vitamin Mix, Teklad (40060)	10.0
Mineral Mix, AIN-76 (170915)	35.0
Calcium Phosphate, dibasic	31.17
Calcium Carbonate	14.05

Footnote

Referred to in the literature as a rescue diet for mice lacking a functional vitamin D receptor. Contains approximately 2200 IU vitamin D/kg diet, along with 20% lactose, 2% calcium, and 1.25% phosphorus. TD.94112 is a version with green color coding. TD.87095 is a version without vitamin D that is colored brown. For questions about this or other vitamin and mineral adjusted diets askanutritionist@inotivco.com.

Selected Nutrient Information ¹

	% by weight	% kcal from
Protein	17.6	19.8
Carbohydrate	59.8	67.5
Fat	5.0	12.7
Kcal/g	3.5	

¹ Values are calculated from ingredient analysis or manufacturer data

Collaborate with a Nutritionist

- askanutritionist@inotivco.com

Teklad diets are designed & manufactured for research purposes only
Teklad™ Custom Diets, Madison, WI | inotivco.com



Key Features

- Purified Diet
- Lactose
- Calcium
- Phosphorus

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotivco.com
- 800.483.5523

Price Quotes | Order Status

- teklad@inotivco.com